

BUFFER TABLETS

convenient method for making buffer solutions
at given concentrations

www.lobachemie.com
info@lobachemie.com



LOBA
Chemie

LABORATORY REAGENTS
& FINE CHEMICALS

Loba Chemie offer a variety of pH buffer tablets for making a buffer solution at a given concentration. These buffer tablets are in complement to our already existing extensive range of buffer line.

- each tablet individually wrapped for a long shelf life
- rapid and convenient method of preparing buffer solution
- perfect for small volumes of pH buffers needed for calibration

Range of buffer tablets

BUFFER TABLETS PH 4.0

(one tablet dissolved in 100 ml of distilled water will give a solution having ph 4.0 at 20°C)

Product Code : 02374

Available packing: 50 Tab, 100 Tab, 1000 Tab

BUFFER TABLETS PH 6.0

(one tablet dissolved in 100 ml of distilled water will give a solution having ph 6.0 at 20°C)

Product Code : 2374D

Available packing: 50 Tab, 100 Tab, 1000 Tab

BUFFER TABLETS PH 6.40

(one tablet dissolved in 100 ml of distilled water will give a solution having ph6.40 at 20°C)

Product Code : 2374G

Available packing: 50 Tab, 100 Tab

BUFFER TABLETS PH 6.40

(one tablet dissolved in 1000 ml solution)
(for preparation of buffer acc. to WEISE staining of blood smear)

Product Code : 2374H

Available packing: 50 Tab, 100 Tab

BUFFER TABLETS PH 6.80

(one tablet dissolved in 100 ml of distilled water will give a solution having ph 6.8 at 20°C)

Product Code : 2374E

Available packing: 50 Tab, 100 Tab, 1000 Tab

BUFFER TABLETS PH 6.80

(one tablet dissolved in 1000 ml solution)
(for preparation of buffer acc. to WEISE staining of blood smear)

Product Code : 2374F

Available packing: 50 Tab, 100 Tab

BUFFER TABLETS PH 7.0

(one tablet dissolved in 100 ml of distilled water will give a solution having ph 7.0 at 20°C)

Product Code : 02375

Available packing: 50 Tab, 100 Tab, 1000 Tab

BUFFER TABLETS PH 7.0

(pH after dissolving 10 tablets in distilled water to make 100 ml solution (at 20°C) : 7.0+/- 0.05)

Product Code : 2375A

Available packing: 50 Tab, 100 Tab

BUFFER TABLETS PH 7.0

(one tablet dissolved in 1000 ml solution)
(for preparation of buffer acc. to WEISE staining of blood smear)

Product Code : 2375B

Available packing: 50 Tab, 100 Tab

BUFFER TABLETS PH 7.20

(one tablet dissolved in 1000 ml solution)
(for preparation of buffer acc. to WEISE staining of blood smear)

Product Code : 2375E

Available packing: 50 Tab, 100 Tab

BUFFER TABLETS PH 9.2

(one tablet dissolved in 100 ml of distilled water will give a solution having ph 9.2 at 20°C)

Product Code : 02376

Available packing: 50 Tab, 100 Tab, 1000 Tab